Easy iMatrix-511

No dilution necessary

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- Ido H et al. The requirement of the glutamic acid residue at the third position from the carboxyl termini of the laminin γ chains in integrin binding by laminins. *Journal of Biological Chemistry*, **282**(15), 11144-11154. 2007
- Taniguchi Y et al. The C-terminal region of laminin β chains modulates the integrin binding affinities of laminins. *Journal of Biological Chemistry*, **284**(12), 7820-7831. 2009
- Miyazaki T et al. Laminin E8 fragments support efficient adhesion and expansion of dissociated human pluripotent stem cells. *Nature communications*, **3**, 1236. 2012
- Nakagawa M et al. A novel efficient feeder-free culture system for the derivation of human induced pluripotent stem cells. *Scientific Reports*, **4**, 3594. 2014
- Doi D et al. Isolation of human induced pluripotent stem cell-derived dopaminergic progenitors by cell sorting for successful transplantation. *Stem cell reports*, **2**(3), 337-350. 2014
- Takahashi Y et al. Resetting transcription factor control circuitry toward ground-state pluripotency in human. *Cell*, **158**(6), 1254-1269. 2014
- Fukuta M et al. Derivation of mesenchymal stromal cells from pluripotent stem cells through a neural crest lineage using small molecule compounds with defined media. *PLoS One*, **9**(12), e112291. 2014
- Burridge PW et al. Chemically defined generation of human cardiomyocytes. *Nature methods*, **11**(8), 855-860. 2014
- Okumura N et al. Laminin-511 and -521 enable efficient in vitro expansion of human corneal endothelial cells laminin-511 and -521 enable expansion of HCECs. *Investigative ophthalmology & visual science*, **56**(5), 2933-2942. 2015
- Sasaki K et al. Robust in vitro induction of human germ cell fate from pluripotent stem cells. *Cell stem cell*, **17**(2), 178-194. 2015
- Hayashi R et al. Co-ordinated ocular development from human iPS cells and recovery of corneal function. *Nature*, **531**(7594), 376-380. 2016
- Matsuno K et al. Redefining definitive endoderm subtypes by robust induction of human induced pluripotent stem cells. *Differentiation*, **92**(5), 281-290. 2016
- Nishimura K et al. Estradiol facilitates functional integration of iPSC-derived dopaminergic neurons into striatal neuronal circuits via activation of integrin $\alpha 5 \beta 1$. *Stem cell reports*, **6**(4), 511-524. 2016
- Takayama K et al. Laminin 411 and 511 promote the cholangiocyte differentiation of human induced pluripotent stem cells. *Biochemical and biophysical research communications*, **474**(1), 91-96. 2016
- Samata B et al. Purification of functional human ES and iPSC-derived midbrain dopaminergic progenitors using LRTM1. *Nature communications*, **7**, 13097. 2016
- Ishikawa T et al. Genetic and pharmacological correction of aberrant dopamine synthesis using patient iPSCs with BH4 metabolism disorders. *Human Molecular Genetics*, ddw339. 2016
- Miyazaki T et al. Efficient adhesion culture of human pluripotent stem cells using laminin fragments in an uncoated manner. *Scientific Reports*, **7**, 41165. 2017
- Goparaju, S. K et al. Rapid differentiation of human pluripotent stem cells into functional neurons by mRNAs encoding transcription factors. *Scientific Reports*, **7**, 42367. 2017

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